

— ABLYBIO, Help Your Research



Phospho-Integrin alpha 4 (Ser1027) (YD12887) Rabbit mAb

货号: **AYD14055**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	115kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat
细胞定位	Membrane
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	The gene encodes a member of the integrin alpha chain family of proteins. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain that function in cell surface adhesion and signaling. The encoded preproprotein is proteolytically processed to generate light and heavy chains that comprise the alpha 4 subunit. This subunit associates with a beta 1 or beta 7 subunit to form an integrin that may play a role in cell motility and migration. This integrin is a therapeutic target for the treatment of multiple sclerosis, Crohn's disease and inflammatory bowel disease. Alternative splicing results in multiple transcript variants.
基因ID	3676
基因名	ITGA4
Swiss	P13612 (https://www.uniprot.org/uniprotkb/P13612/entry)
别名	Phospho-Integrin alpha 4 (Ser1027) (YD12887), Phospho-Integrin alpha 4 (Ser1027) (YD12887) Rabbit mAb, ITGA4, CD49 antigen-like family member D, Integrin alpha-IV, VLA-4 subunit alpha, CD49D

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)