

— ABLYBIO, Help Your Research



Dynein Light Chain Tctex Type 1 (YD32413) Rabbit mAb

货号: **AYD14065**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	12kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,MCF7,Mouse thymus,Mouse heart
细胞定位	Golgi apparatus, Cytoplasm, cytoskeleton, spindle
纯化	亲和纯化


抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a component of the motor complex, cytoplasmic dynein, which transports cellular cargo along microtubules in the cell. The encoded protein regulates the length of primary cilia which are sensory organelles found on the surface of cells. The protein encoded by this gene interacts with viral proteins, like the minor capsid protein L2 of human papillomavirus, and is required for dynein-mediated delivery of the viral nucleic acid to the host nucleus. This protein interacts with oncogenic nucleoporins to disrupt gene regulation and cause leukemic transformation. Pseudogenes of this gene are present on chromosomes 4 and 17. Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	6993
基因名	DYNLT1
Swiss	P63172 (https://www.uniprot.org/uniprotkb/P63172/entry)
别名	Dynein Light Chain Tctex Type 1 (YD32413), Dynein Light Chain Tctex Type 1 (YD32413) Rabbit mAb, DYNLT1, Protein CW-1, T-complex testis-specific protein 1 homolog, TCTEL1, TCTEX-1, TCTEX1

产品验证





实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)