

— ABLYBIO, Help Your Research



Phospho-HSF1 (Ser320) (YD12830) Rabbit mAb

货号: **AYD14067**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	57kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,K-562,HepG2,MCF7
细胞定位	Nucleus, Cytoplasm, nucleoplasm, perinuclear region, cytoskeleton, spindle pole, microtubule organizing center, centrosome, Chromosome, centromere, kinetochores
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The product of this gene is a heat-shock transcription factor. Transcription of heat-shock genes is rapidly induced after temperature stress. Hsp90, by itself and/or associated with multichaperone complexes, is a major repressor of this gene.
------	--

基因ID	3297
基因名	HSF1
Swiss	Q00613 (https://www.uniprot.org/uniprotkb/Q00613/entry)
别名	Phospho-HSF1 (Ser320) (YD12830),Phospho-HSF1 (Ser320) (YD12830) Rabbit mAb,HSF1,Heat shock transcription factor 1,HSTF1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)