

— ABLYBIO, Help Your Research



GSK3 (alpha/beta) (YD14064) Rabbit mAb

货号: **AYD14079**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	51kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-431,HeLa,B-cell,Mouse liver
细胞定位	
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a multifunctional Ser/Thr protein kinase that is implicated in the control of several regulatory proteins including glycogen synthase, and transcription factors, such as JUN. It also plays a role in the WNT and PI3K signaling pathways, as well as regulates the production of beta-amyloid peptides associated with Alzheimer's disease.
------	---

基因ID	2931
基因名	GSK3A
Swiss	P49840 (https://www.uniprot.org/uniprotkb/P49840/entry)
别名	GSK3 (alpha/beta) (YD14064),GSK3 (alpha/beta) (YD14064) Rabbit mAb,GSK3A,Serine/threonine-protein kinase GSK3A

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)