

— ABLYBIO, Help Your Research



# Nucleoside triphosphate phosphohydrolase (YD20100) Rabbit mAb

货号: AYD14095

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC IP
推荐浓度	
理论分子量	21kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,A-431,293T
细胞定位	
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a non-specific nucleoside triphosphatase that is slow-acting in vitro. This gene is overexpressed in many tumor tissues, and while it is not essential for the cell, overexpression is cytotoxic. However, the cytotoxicity is not related to its triphosphatase activity.
基因ID	84284
基因名	NTPCR
Swiss	Q9BSD7 ( <a href="https://www.uniprot.org/uniprotkb/Q9BSD7/entry">https://www.uniprot.org/uniprotkb/Q9BSD7/entry</a> )
别名	Nucleoside triphosphate phosphohydrolase (YD20100),Nucleoside triphosphate phosphohydrolase (YD20100) Rabbit mAb,NTPCR,Nucleoside triphosphate phosphohydrolase,C1orf57

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)