

— ABLYBIO, Help Your Research



Phospho-alpha Synuclein (Ser129) (YD11937) Rabbit mAb

货号: **AYD14104**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	14kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse brain,Rat brain
细胞定位	Cytoplasm, Membrane, Nucleus, Synapse, Secreted, Cell projection, axon
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Alpha-synuclein is a member of the synuclein family, which also includes beta- and gamma-synuclein. Synucleins are abundantly expressed in the brain and alpha- and beta-synuclein inhibit phospholipase D2 selectively. SNCA may serve to integrate presynaptic signaling and membrane trafficking. Defects in SNCA have been implicated in the pathogenesis of Parkinson disease. SNCA peptides are a major component of amyloid plaques in the brains of patients with Alzheimer's disease. Alternatively spliced transcripts encoding different isoforms have been identified for this gene.
基因ID	6622
基因名	SNCA
Swiss	P37840 (https://www.uniprot.org/uniprotkb/P37840/entry)
别名	Phospho-alpha Synuclein (Ser129) (YD11937), Phospho-alpha Synuclein (Ser129) (YD11937) Rabbit mAb, SNCA, Non-A beta component of AD amyloid, Non-A4 component of amyloid precursor, NACP, PARK1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)