

— ABLYBIO, Help Your Research



Phospho-Lyn (Tyr507) (YD11595) Rabbit mAb

货号: **AYD14116**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	59kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HeLa,K-562,Mouse lung,Rat lung
细胞定位	Cell membrane, Nucleus, Cytoplasm, perinuclear region, Golgi apparatus, Membrane
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a tyrosine protein kinase, which maybe involved in the regulation of mast cell degranulation, and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	4067
基因名	LYN
Swiss	P07948 (https://www.uniprot.org/uniprotkb/P07948/entry)
别名	Phospho-Lyn (Tyr507) (YD11595), Phospho-Lyn (Tyr507) (YD11595) Rabbit mAb, LYN, Lck/Yes-related novel protein tyrosine kinase, V-yes-1 Yamaguchi sarcoma viral related oncogene homolog, p53Lyn, p56Lyn, JTK8

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)