

— ABLYBIO, Help Your Research



# Phospho-EGFR (Ser1070/Ser1071) (YD11572) Rabbit mAb

货号: AYD14120

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF IP
推荐浓度	
理论分子量	134kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cell membrane, Endoplasmic reticulum membrane, Golgi apparatus membrane, Nucleus membrane, Endosome, Endosome membrane, Nucleus, Secreted
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a transmembrane glycoprotein that is a member of the protein kinase superfamily. This protein is a receptor for members of the epidermal growth factor family. EGFR is a cell surface protein that binds to epidermal growth factor. Binding of the protein to a ligand induces receptor dimerization and tyrosine autophosphorylation and leads to cell proliferation. Mutations in this gene are associated with lung cancer.
基因ID	1956
基因名	EGFR
Swiss	P00533 ( <a href="https://www.uniprot.org/uniprotkb/P00533/entry">https://www.uniprot.org/uniprotkb/P00533/entry</a> )
别名	Phospho-EGFR (Ser1070/Ser1071) (YD11572),Phospho-EGFR (Ser1070/Ser1071) (YD11572) Rabbit mAb,EGFR,Proto-oncogene c-ErbB-1,Receptor tyrosine-protein kinase erbB-1,ERBB,ERBB1,HER1

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)