

— ABLYBIO, Help Your Research



SORL1 (YD21036) Rabbit mAb

货号: **AYD14133**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	248kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y
细胞定位	Golgi apparatus membrane, Golgi apparatus, trans-Golgi network membrane, Endosome membrane, Early endosome membrane, Recycling endosome membrane, Endoplasmic reticulum membrane, Endosome, multivesicular body membrane, Cell membrane, Cytoplasmic vesicle, secretory vesicle membrane, Secreted
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a mosaic protein that belongs to at least two families: the vacuolar protein sorting 10 (VPS10) domain-containing receptor family, and the low density lipoprotein receptor (LDLR) family. The encoded protein also contains fibronectin type III repeats and an epidermal growth factor repeat. The encoded preproprotein is proteolytically processed to generate the mature receptor, which likely plays roles in endocytosis and sorting. Mutations in this gene may be associated with Alzheimer's disease.
基因ID	6653
基因名	SORL1
Swiss	Q92673 (https://www.uniprot.org/uniprotkb/Q92673/entry)
别名	SORL1 (YD21036),SORL1 (YD21036) Rabbit mAb,SORL1,Low-density lipoprotein receptor relative with 11 ligand-binding repeats,SorLA-1,Sorting protein-related receptor containing LDLR class A repeats,C11orf32

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)