

— ABLYBIO, Help Your Research



# EpCAM (YD20941) Rabbit mAb

货号: **AYD14173**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF IP
推荐浓度	
理论分子量	35kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Lateral cell membrane, Cell junction, tight junction
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	May act as a physical homophilic interaction molecule between intestinal epithelial cells (IECs) and intraepithelial lymphocytes (IELs) at the mucosal epithelium for providing immunological barrier as a first line of defense against mucosal infection. Plays a role in embryonic stem cells proliferation and differentiation. Up-regulates the expression of FABP5, MYC and cyclins A and E (By similarity)
------	---

基因ID	3134
基因名	Epcam
Swiss	Q99JW5 ( <a href="https://www.uniprot.org/uniprotkb/Q99JW5/entry">https://www.uniprot.org/uniprotkb/Q99JW5/entry</a> )
别名	EpCAM (YD20941),EpCAM (YD20941) Rabbit mAb,Epcam,Epithelial glycoprotein 314,Protein 289A,Tumor-associated calcium signal transducer 1,Tacstd1

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)