

— ABLYBIO, Help Your Research



TGF beta Receptor II (YD20933) Rabbit mAb

货号: **AYD14182**

产品信息

| | |
|-------|---|
| 反应 | Human |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | |
| 理论分子量 | 65kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Mouse uterus |
| 细胞定位 | Cell membrane, Membrane raft, Secreted |
| 纯化 | 亲和纯化 |


抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|--|
| 研究背景 | This gene encodes a member of the Ser/Thr protein kinase family and the TGFβ receptor subfamily. The encoded protein is a transmembrane protein that has a protein kinase domain, forms a heterodimeric complex with another receptor protein, and binds TGF-β. This receptor/ligand complex phosphorylates proteins, which then enter the nucleus and regulate the transcription of a subset of genes related to cell proliferation. Mutations in this gene have been associated with Marfan Syndrome, Loeys-Deitz Aortic Aneurysm Syndrome, and the development of various types of tumors. Alternatively spliced transcript variants encoding different isoforms have been characterized. |
| 基因ID | 7048 |
| 基因名 | TGFBR2 |
| Swiss | P37173 (https://www.uniprot.org/uniprotkb/P37173/entry) |
| 别名 | TGF beta Receptor II (YD20933),TGF beta Receptor II (YD20933) Rabbit mAb,TGFBR2,TGF-beta type II receptor,Transforming growth factor-beta receptor type II |

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)