

— ABLYBIO, Help Your Research



# FGL1 (YD20886) Rabbit mAb

货号: **AYD14209**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	36kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	Secreted
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	Fibrinogen-like protein 1 (FGL1), a liver-secreted protein, is a major LAG-3 functional ligand independent from MHC-II. The FGL1-LAG-3 pathway maybe an important immune evasion mechanism and could contribute to current cancer immunotherapy efforts for several reasons. FGL1 inhibits antigen-specific T cell activation, and ablation of FGL1 in mice promotes T cell immunity. Blockade of the FGL1-LAG-3 interaction by monoclonal antibodies stimulates tumor immunity and is therapeutic against established mouse tumors in a receptor-ligand inter-dependent manner. FGL1 is highly produced by human cancer cells, and elevated FGL1 in the plasma of cancer patients is associated with a poor prognosis and resistance to anti-PD-1/B7-H1 therapy
基因ID	2267
基因名	FGL1
Swiss	Q08830 ( <a href="https://www.uniprot.org/uniprotkb/Q08830/entry">https://www.uniprot.org/uniprotkb/Q08830/entry</a> )
别名	FGL1 (YD20886),FGL1 (YD20886) Rabbit mAb,FGL1,HP-041,Hepassocin,Hepatocyte-derived fibrinogen-related protein 1,Liver fibrinogen-related protein 1,HFREP1

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)