

— ABLYBIO, Help Your Research



IL-10 (YD20530) Rabbit mAb

货号: **AYD14268**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	21kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	RAW264.7
细胞定位	Secreted
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a cytokine produced primarily by monocytes and to a lesser extent by lymphocytes. This cytokine has pleiotropic effects in immunoregulation and inflammation. It down-regulates the expression of Th1 cytokines, MHC class II Ags, and costimulatory molecules on macrophages. It also enhances B cell survival, proliferation, and antibody production. This cytokine can block NF-kappa B activity, and is involved in the regulation of the JAK-STAT signaling pathway. Knockout studies in mice suggested the function of this cytokine as an essential immunoregulator in the intestinal tract. Mutations in this gene are associated with an increased susceptibility to HIV-1 infection and rheumatoid arthritis.
基因ID	3586
基因名	IL10
Swiss	P22301 (https://www.uniprot.org/uniprotkb/P22301/entry)
别名	IL-10 (YD20530),IL-10 (YD20530) Rabbit mAb,IL10,Cytokine synthesis inhibitory factor

产品验证

RAW264.7+5ul LPS



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)