

— ABLYBIO, Help Your Research



Fos B (YD20466) Rabbit mAb

货号: **AYD14279**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF IP
推荐浓度	
理论分子量	36kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse brain,Rat brain
细胞定位	Nucleus
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	2354
基因名	FOSB
Swiss	P53539 (https://www.uniprot.org/uniprotkb/P53539/entry)
别名	Fos B (YD20466),Fos B (YD20466) Rabbit mAb,FOSB,FosB proto-oncogene,AP-1 transcription factor subunit,G0/G1 switch regulatory protein 3,Transcription factor AP-1 subunit FosB,G0S3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)