

— ABLYBIO, Help Your Research



# SERPINB1 (YD20438) Rabbit mAb

货号: **AYD14284**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	43kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,K-562,Mouse lung,Rat lung,Rat spleen
细胞定位	Secreted, Cytoplasm, Cytolytic granule, Early endosome
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the serpin family of proteinase inhibitors. Members of this family maintain homeostasis by neutralizing overexpressed proteinase activity through their function as suicide substrates. This protein inhibits the neutrophil-derived proteinases neutrophil elastase, cathepsin G, and proteinase-3 and thus protects tissues from damage at inflammatory sites. Alternative splicing results in multiple transcript variants.
基因ID	1992
基因名	SERPINB1
Swiss	P30740 ( <a href="https://www.uniprot.org/uniprotkb/P30740/entry">https://www.uniprot.org/uniprotkb/P30740/entry</a> )
别名	SERPINB1 (YD20438),SERPINB1 (YD20438) Rabbit mAb,SERPINB1,Monocyte/neutrophil elastase inhibitor, Peptidase inhibitor 2,Serpin B1,ELANH2,MNEI,PI2

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn//www.ablybio.cn>)