

— ABLYBIO, Help Your Research



RBMX2 (YD20326) Rabbit mAb

货号: **AYD14299**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	37kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus
纯化	亲和纯化

抗原信息



抗原信息	
------	--

靶点信息

研究背景	Involved in pre-mRNA splicing as component of the activated spliceosome. As a component of the minor spliceosome, involved in the splicing of U12-type introns in pre-mRNAs (Probable)
基因ID	51634

基因名	RBMX2
Swiss	Q9Y388 (https://www.uniprot.org/uniprotkb/Q9Y388/entry)
别名	RBMX2 (YD20326),RBMX2 (YD20326) Rabbit mAb,RBMX2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)