

— ABLYBIO, Help Your Research



# GGA2 (YD20262) Rabbit mAb

货号: **AYD14311**

## 产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	67kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Raji,MCF7,THP-1,22Rv1,Mouse brain,Mouse kidney
细胞定位	Golgi apparatus, trans-Golgi network membrane, Endosome membrane, Early endosome membrane
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the Golgi-localized, gamma adaptin ear-containing, ARF-binding (GGA) family. This family includes ubiquitous coat proteins that regulate the trafficking of proteins between the trans-Golgi network and the lysosome. These proteins share an amino-terminal VHS domain which mediates sorting of the mannose 6-phosphate receptors at the trans-Golgi network. They also contain a carboxy-terminal region with homology to the ear domain of gamma-adaptins. This family member may play a significant role in cargo molecules regulation and clathrin-coated vesicle assembly.
基因ID	23062
基因名	GGA2
Swiss	Q9UJY4 ( <a href="https://www.uniprot.org/uniprotkb/Q9UJY4/entry">https://www.uniprot.org/uniprotkb/Q9UJY4/entry</a> )
别名	GGA2 (YD20262),GGA2 (YD20262) Rabbit mAb,GGA2,Gamma-adaptin-related protein 2,Golgi-localized,gamma ear-containing,ARF-binding protein 2,VHS domain and ear domain of gamma-adaptin,KIAA1080

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)