

— ABLYBIO, Help Your Research



SEC24D (YD31537) Rabbit mAb

货号: **AYD14312**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	113kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,BxPC-3,A-375,DU145,Mouse liver
细胞定位	Cytoplasmic vesicle, COPII-coated vesicle membrane, Endoplasmic reticulum membrane, Cytoplasm, cytosol
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a member of the SEC24 subfamily of the SEC23/SEC24 family, which is involved in vesicle trafficking. The encoded protein has similarity to yeast Sec24p component of COPII. COPII is the coat protein complex responsible for vesicle budding from the ER. This gene product is implicated in the shaping of the vesicle, and also in cargo selection and concentration. Mutations in this gene have been associated with Cole-Carpenter syndrome, a disorder affecting bone formation, resulting in craniofacial malformations and bones that break easily. Alternative splicing results in multiple transcripts encoding different isoforms.
基因ID	9871
基因名	SEC24D
Swiss	O94855 (https://www.uniprot.org/uniprotkb/O94855/entry)
别名	SEC24D (YD31537),SEC24D (YD31537) Rabbit mAb,SEC24D,SEC24-related protein D,KIAA0755

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)