

— ABLYBIO, Help Your Research



MIP-3 beta (YD20242) Rabbit mAb

货号: **AYD14322**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	11kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Secreted
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This antimicrobial gene is one of several CC cytokine genes clustered on the p-arm of chromosome 9. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterized by two adjacent cysteines. The cytokine encoded by this gene may play a role in normal lymphocyte recirculation and homing. It also plays an important role in trafficking of T cells in thymus, and in T cell and B cell migration to secondary lymphoid organs. It specifically binds to chemokine receptor CCR7.
基因ID	6363
基因名	CCL19
Swiss	Q99731 (https://www.uniprot.org/uniprotkb/Q99731/entry)
别名	MIP-3 beta (YD20242),MIP-3 beta (YD20242) Rabbit mAb,CCL19,Beta-chemokine exodus-3,CK beta-11,Epstein-Barr virus-induced molecule 1 ligand chemokine,Macrophage inflammatory protein 3 beta,Small-inducible cytokine A19,ELC,MIP3B,SCYA19

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)