

— ABLYBIO, Help Your Research



# IKK alpha/beta (YD20221) Rabbit mAb

货号: **AYD14327**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	87kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse spleen
细胞定位	Cytoplasm, Nucleus, Membrane raft
纯化	亲和纯化

## 抗原信息



抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene phosphorylates the inhibitor in the inhibitor/NF-kappa-B complex, causing dissociation of the inhibitor and activation of NF-kappa-B. The encoded protein itself is found in a complex of proteins. Several transcript variants, some protein-coding and some not, have been found for this gene.
------	--

基因ID	3551
基因名	IKBKB
Swiss	O14920 ( <a href="https://www.uniprot.org/uniprotkb/O14920/entry">https://www.uniprot.org/uniprotkb/O14920/entry</a> )
别名	IKK alpha/beta (YD20221),IKK alpha/beta (YD20221) Rabbit mAb,IKBKB,I-kappa-B kinase 2,Nuclear factor NF-kappa-B inhibitor kinase beta,Serine/threonine protein kinase IKBKB,IKKB

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)