

— ABLYBIO, Help Your Research



# CD11b (YD20208) Rabbit mAb

货号: **AYD14335**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF IP
推荐浓度	
理论分子量	127kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	THP-1,U-937,Mouse lung,Mouse brain,Mouse spleen,Rat lung
细胞定位	Cell membrane, Membrane raft
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes the integrin alpha M chain. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This I-domain containing alpha integrin combines with the beta 2 chain (ITGB2) to form a leukocyte-specific integrin referred to as macrophage receptor 1 ('Mac-1'), or inactivated-C3b (iC3b) receptor 3 ('CR3'). The alpha M beta 2 integrin is important in the adherence of neutrophils and monocytes to stimulated endothelium, and also in the phagocytosis of complement coated particles. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	3684
基因名	ITGAM
Swiss	P11215 ( <a href="https://www.uniprot.org/uniprotkb/P11215/entry">https://www.uniprot.org/uniprotkb/P11215/entry</a> )
别名	CD11b (YD20208),CD11b (YD20208) Rabbit mAb,ITGAM,CD11 antigen-like family member B,CR-3 alpha chain,Cell surface glycoprotein MAC-1 subunit alpha,Leukocyte adhesion receptor MO1,Neutrophil adherence receptor,CD11B,CR3A

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)