

— ABLYBIO, Help Your Research



# NCX1 (YD20194) Rabbit mAb

货号: **AYD14342**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	109kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,RD,SH-SY5Y,Mouse brain
细胞定位	Cell membrane
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	In cardiac myocytes, Ca(2+) concentrations alternate between high levels during contraction and low levels during relaxation. The increase in Ca(2+) concentration during contraction is primarily due to release of Ca(2+) from intracellular stores. However, some Ca(2+) also enters the cell through the sarcolemma (plasma membrane). During relaxation, Ca(2+) is sequestered within the intracellular stores. To prevent overloading of intracellular stores, the Ca(2+) that entered across the sarcolemma must be extruded from the cell. The Na(+)-Ca(2+) exchanger is the primary mechanism by which the Ca(2+) is extruded from the cell during relaxation. In the heart, the exchanger may play a key role in digitalis action. The exchanger is the dominant mechanism in returning the cardiac myocyte to its resting state following excitation.
基因ID	6546
基因名	SLC8A1
Swiss	P32418 ( <a href="https://www.uniprot.org/uniprotkb/P32418/entry">https://www.uniprot.org/uniprotkb/P32418/entry</a> )
别名	NCX1 (YD20194),NCX1 (YD20194) Rabbit mAb,SLC8A1,Na(+)/Ca(2+)-exchange protein 1,Solute carrier family 8 member 1,CNC,NCX1

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)