

— ABLYBIO, Help Your Research



BMI1 (YD19989) Rabbit mAb

货号: **AYD14345**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	37kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,K-562,Mouse heart,Mouse thymus,Rat thymus,Rat testis
细胞定位	Nucleus, Cytoplasm
纯化	亲和纯化




抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a ring finger protein that is major component of the polycomb group complex 1 (PRC 1). This complex functions through chromatin remodeling as an essential epigenetic repressor of multiple regulatory genes involved in embryonic development and self-renewal in somatic stem cells. This protein also plays a central role in DNA damage repair. This gene is an oncogene and aberrant expression is associated with numerous cancers and is associated with resistance to certain chemotherapies. A pseudogene of this gene is found on chromosome X. Read-through transcription also exists between this gene and the upstream COMM domain containing 3 (COMMD3) gene.
基因ID	100532731
基因名	BMI1
Swiss	P35226 (https://www.uniprot.org/uniprotkb/P35226/entry)
别名	BMI1 (YD19989),BMI1 (YD19989) Rabbit mAb,BMI1,Polycomb group RING finger protein 4,RING finger protein 51,PCGF4,RNF51

产品验证





实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)