

— ABLYBIO, Help Your Research



Oct-2 (YD33484) Rabbit mAb

货号: **AYD14389**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	51kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Raji, Daudi, HepG2, NIH/3T3
细胞定位	Cytoplasm, Nucleus
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a homeobox-containing transcription factor of the POU domain family. The encoded protein binds the octamer sequence 5'-ATTGTCAT-3', a common transcription factor binding site in immunoglobulin gene promoters. Several transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	5452
基因名	POU2F2
Swiss	P09086 (https://www.uniprot.org/uniprotkb/P09086/entry)
别名	Oct-2 (YD33484),Oct-2 (YD33484) Rabbit mAb,POU2F2,Lymphoid-restricted immunoglobulin octamer-binding protein NF-A2,Octamer-binding protein 2,Octamer-binding transcription factor 2,OCT2,OTF2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)