

— ABLYBIO, Help Your Research



Pirh2 (YD11887) Rabbit mAb

货号: **AYD14422**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF IP
推荐浓度	
理论分子量	30kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HCT116
细胞定位	Nucleus, Nucleus speckle, Cytoplasm
纯化	亲和纯化

抗原信息



抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene has ubiquitin ligase activity. It mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including tumor protein 53, histone deacetylase 1, and cyclin-dependent kinase inhibitor 1B, thus regulating their levels and cell cycle progression. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.
------	--

基因ID	25898
基因名	RCHY1
Swiss	Q96PM5 (https://www.uniprot.org/uniprotkb/Q96PM5/entry)
别名	Pirh2 (YD11887), Pirh2 (YD11887) Rabbit mAb, RCHY1, Androgen receptor N-terminal-interacting protein, C H-rich-interacting match with PLAG1, E3 ubiquitin-protein ligase Pirh2, RING finger protein 199, RING-type E3 ubiquitin transferase RCHY1, Zinc finger protein 363

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)