

RASSF2 (YD12747) Rabbit mAb

货号: **AYD14442**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	38kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Rat lung
细胞定位	Nucleus, Cytoplasm, Chromosome, centromere, kinetochores
纯化	


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a protein that contains a Ras association domain. Similar to its cattle and sheep counterparts, this gene is located near the prion gene. Two alternatively spliced transcripts encoding the same isoform have been reported.
基因ID	9770
基因名	RASSF2
Swiss	P50749
别名	RASSF2 (YD12747)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn