

Vitronectin (YD19605) Rabbit mAb

货号: **AYD14487**

产品信息

| | |
|-------|---|
| 反应 | Human |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC-P |
| 推荐浓度 | |
| 理论分子量 | 54kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Human serum,Rat serum |
| 细胞定位 | Secreted, extracellular space, Parasitophorous vacuole |
| 纯化 | |

抗原信息


| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|--|
| 研究背景 | The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond. |
| 基因ID | 7448 |
| 基因名 | VTN |
| Swiss | P04004 |

| | |
|----|-----------------------|
| 别名 | Vitronectin (YD19605) |
|----|-----------------------|

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn