

FOXF1 (YD19581) Rabbit mAb

货号: **AYD14488**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	40kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse stomach
细胞定位	Nucleus
纯化	



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in the regulation of pulmonary genes as well as embryonic development.
基因ID	2294
基因名	FOXF1
Swiss	Q12946
别名	FOXF1 (YD19581)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn