

FAK (YD32081) Rabbit mAb

货号: **AYD14504**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	119kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	OVCAR3,Mouse testis,Mouse brain,Rat testis
细胞定位	Cell junction, focal adhesion, Cell membrane, Cytoplasm, perinuclear region, cell cortex, cytoskeleton, microtubule organizing center, centrosome, Nucleus, cilium basal body
纯化	

抗原信息



抗原信息	
------	--

靶点信息

研究背景	This gene encodes a cytoplasmic protein tyrosine kinase which is found concentrated in the focal adhesions that form between cells growing in the presence of extracellular matrix constituents. The encoded protein is a member of the FAK subfamily of protein tyrosine kinases but lacks significant sequence similarity to kinases from other subfamilies. Activation of this gene may be an important early step in cell growth and intracellular signal transduction pathways triggered in response to certain neural peptides or to cell interactions with the extracellular matrix. Several transcript variants encoding different isoforms have been found for this gene.
基因ID	5747

基因名	PTK2
Swiss	Q05397
别名	FAK (YD32081)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn