

SIRT2 (YD19368) Rabbit mAb

货号: **AYD14514**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC IP
推荐浓度	
理论分子量	43kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Nucleus, Cytoplasm, perinuclear region, cytoskeleton, microtubule organizing center, centrosome, centriole, spindle, Midbody, Chromosome, Perikaryon, Cell projection, growth cone, Myelin membrane
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the sirtuin family of proteins, homologs to the yeast Sir2 protein. Members of the sirtuin family are characterized by a sirtuin core domain and grouped into four classes. The functions of human sirtuins have not yet been determined; however, yeast sirtuin proteins are known to regulate epigenetic gene silencing and suppress recombination of rDNA. Studies suggest that the human sirtuins may function as intracellular regulatory proteins with mono-ADP-ribosyltransferase activity. The protein encoded by this gene is included in class I of the sirtuin family. Several transcript variants are resulted from alternative splicing of this gene.
基因ID	22933

基因名	SIRT2
Swiss	Q8IXJ6
别名	SIRT2 (YD19368)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn