

— ABLYBIO, Help Your Research



# ACTN3 (YD19571) Rabbit mAb

货号: **AYD14518**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC
推荐浓度	
理论分子量	103kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,Mouse skeletal muscle,Rat skeletal muscle
细胞定位	
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the alpha-actin binding protein gene family. The encoded protein is primarily expressed in skeletal muscle and functions as a structural component of sarcomeric Z line. This protein is involved in crosslinking actin containing thin filaments. An allelic polymorphism in this gene results in both coding and non-coding variants; the reference genome represents the coding allele. The non-functional allele of this gene is associated with elite athlete status.
基因ID	89
基因名	ACTN3
Swiss	Q08043 ( <a href="https://www.uniprot.org/uniprotkb/Q08043/entry">https://www.uniprot.org/uniprotkb/Q08043/entry</a> )
别名	ACTN3 (YD19571),ACTN3 (YD19571) Rabbit mAb,ACTN3,Alpha-actinin skeletal muscle isoform 3,F-actin cross-linking protein

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)