

— ABLYBIO, Help Your Research



# Netrin 1 (YD19045) Rabbit mAb

货号: **AYD14569**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	68kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,HepG2
细胞定位	Secreted, Cytoplasm
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	Netrin is included in a family of laminin-related secreted proteins. The function of this gene has not yet been defined; however, netrin is thought to be involved in axon guidance and cell migration during development. Mutations and loss of expression of netrin suggest that variation in netrin may be involved in cancer development.
------	---

基因ID	9423
基因名	NTN1
Swiss	O95631 ( <a href="https://www.uniprot.org/uniprotkb/O95631/entry">https://www.uniprot.org/uniprotkb/O95631/entry</a> )
别名	Netrin 1 (YD19045),Netrin 1 (YD19045) Rabbit mAb,NTN1,Epididymis tissue protein Li 131P,NTN1L

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)