

# FAM50A (YD14853) Rabbit mAb

货号: **AYD14590**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	40kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,PC-12,NIH/3T3,Mouse placenta
细胞定位	Nucleus
纯化	





## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene belongs to the FAM50 family. The encoded protein is highly conserved in length and sequence across different species. It is a basic protein containing a nuclear localization signal, and may function as a DNA-binding protein or a transcriptional factor. [provided by RefSeq, Sep 2009]
基因ID	9130
基因名	FAM50A
Swiss	Q14320
别名	FAM50A (YD14853)

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)