

— ABLYBIO, Help Your Research



# PRPF8 (YD20286) Rabbit mAb

货号: **AYD14628**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	274kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat
细胞定位	Nucleus, Nucleus speckle
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	Pre-mRNA splicing occurs in 2 sequential transesterification steps. The protein encoded by this gene is a component of both U2- and U12-dependent spliceosomes, and found to be essential for the catalytic step II in pre-mRNA splicing process. It contains several WD repeats, which function in protein-protein interactions. This protein has a sequence similarity to yeast Prp8 protein. This gene is a candidate gene for autosomal dominant retinitis pigmentosa.
基因ID	10594
基因名	PRPF8
Swiss	Q6P2Q9 ( <a href="https://www.uniprot.org/uniprotkb/Q6P2Q9/entry">https://www.uniprot.org/uniprotkb/Q6P2Q9/entry</a> )
别名	PRPF8 (YD20286),PRPF8 (YD20286) Rabbit mAb,PRPF8,220 kDa U5 snRNP-specific protein,PRP8 homolog, Splicing factor Prp8,p220,PRPC8

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)