

NSDHL (YD32821) Rabbit mAb

货号: **AYD14630**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	42kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Endoplasmic reticulum membrane, Lipid droplet
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is localized in the endoplasmic reticulum and is involved in cholesterol biosynthesis. Mutations in this gene are associated with CHILD syndrome, which is a X-linked dominant disorder of lipid metabolism with disturbed cholesterol biosynthesis, and typically lethal in males. Alternatively spliced transcript variants with differing 5' UTR have been found for this gene.
基因ID	50814
基因名	NSDHL
Swiss	Q15738

别名	NSDHL (YD32821)
----	-----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn