

— ABLYBIO, Help Your Research



# IL-6 (YD18791) Rabbit mAb

货号: **AYD14638**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF IP
推荐浓度	
理论分子量	24kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Secreted
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a cytokine that functions in inflammation and the maturation of B cells. In addition, the encoded protein has been shown to be an endogenous pyrogen capable of inducing fever in people with autoimmune diseases or infections. The protein is primarily produced at sites of acute and chronic inflammation, where it is secreted into the serum and induces a transcriptional inflammatory response through interleukin 6 receptor, alpha. The functioning of this gene is implicated in a wide variety of inflammation-associated disease states, including susceptibility to diabetes mellitus and systemic juvenile rheumatoid arthritis. Alternative splicing results in multiple transcript variants.
基因ID	3569
基因名	IL6
Swiss	P05231 ( <a href="https://www.uniprot.org/uniprotkb/P05231/entry">https://www.uniprot.org/uniprotkb/P05231/entry</a> )
别名	IL-6 (YD18791),IL-6 (YD18791) Rabbit mAb,IL6,B-cell stimulatory factor 2,CTL differentiation factor,Hybridoma growth factor,Interferon beta-2,IFNB2

## 产品验证

Recombinant Human IL-6 protein



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)