

— ABLYBIO, Help Your Research



# Semaphorin 7A (YD20676) Rabbit mAb

货号: **AYD14643**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC IP
推荐浓度	
理论分子量	75kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HT-29,Raji
细胞定位	Cell membrane
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the semaphorin family of proteins. The encoded preproprotein is proteolytically processed to generate the mature glycosylphosphatidylinositol (GPI)-anchored membrane glycoprotein. The encoded protein is found on activated lymphocytes and erythrocytes and may be involved in immunomodulatory and neuronal processes. The encoded protein carries the John Milton Hagen (JMH) blood group antigens. Mutations in this gene may be associated with reduced bone mineral density (BMD). An alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed.
基因ID	8482
基因名	SEMA7A
Swiss	O75326 ( <a href="https://www.uniprot.org/uniprotkb/O75326/entry">https://www.uniprot.org/uniprotkb/O75326/entry</a> )
别名	Semaphorin 7A (YD20676),Semaphorin 7A (YD20676) Rabbit mAb,SEMA7A,CDw108,JMH blood group antigen,John-Milton-Hagen human blood group Ag,Semaphorin-K1,Semaphorin-L,CD108,SEMAL

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)