

— ABLYBIO, Help Your Research



SP1 (YD20064) Rabbit mAb

货号: **AYD14645**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	81kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	Nucleus, Cytoplasm
纯化	亲和纯化




抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a zinc finger transcription factor that binds to GC-rich motifs of many promoters. The encoded protein is involved in many cellular processes, including cell differentiation, cell growth, apoptosis, immune responses, response to DNA damage, and chromatin remodeling. Post-translational modifications such as phosphorylation, acetylation, glycosylation, and proteolytic processing significantly affect the activity of this protein, which can be an activator or a repressor. Three transcript variants encoding different isoforms have been found for this gene.
基因ID	6667
基因名	SP1
Swiss	P08047 (https://www.uniprot.org/uniprotkb/P08047/entry)
别名	SP1 (YD20064),SP1 (YD20064) Rabbit mAb,SP1,TSFP1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)