

— ABLYBIO, Help Your Research



# GRP78 BiP (YD13981) Rabbit mAb

货号: **AYD14659**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	72kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,U-87MG,Mouse liver,Mouse brain,Rat liver
细胞定位	Endoplasmic reticulum lumen, Melanosome, Cytoplasm, Cell surface
纯化	亲和纯化

## 抗原信息



抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the heat shock protein 70 (HSP70) family. It is localized in the lumen of the endoplasmic reticulum (ER), and is involved in the folding and assembly of proteins in the ER. As this protein interacts with many ER proteins, it may play a key role in monitoring protein transport through the cell.
------	--

基因ID	3309
基因名	HSPA5
Swiss	P11021 ( <a href="https://www.uniprot.org/uniprotkb/P11021/entry">https://www.uniprot.org/uniprotkb/P11021/entry</a> )
别名	GRP78 BiP (YD13981),GRP78 BiP (YD13981) Rabbit mAb,HSPA5,78 kDa glucose-regulated protein,Binding-immunoglobulin protein,Heat shock protein 70 family protein 5,Heat shock protein family A member 5,Immunoglobulin heavy chain-binding protein,GRP78

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)