

RanBP3 (YD18561) Rabbit mAb

货号: **AYD14671**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	60kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,LO2,HT-29,SW620,U-251MG
细胞定位	Cytoplasm, Nucleus
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a protein with a RanBD1 domain that is found in both the nucleus and cytoplasm. This protein plays a role in nuclear export as part of a heteromeric complex. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
基因ID	8498
基因名	RANBP3
Swiss	Q9H6Z4
别名	RanBP3 (YD18561)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn