

— ABLYBIO, Help Your Research



EpCAM (YD18511) Rabbit mAb

货号: **AYD14677**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	35kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,Mouse colon,Rat colon
细胞定位	Lateral cell membrane, Cell junction, tight junction
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a carcinoma-associated antigen and is a member of a family that includes at least two type I membrane proteins. This antigen is expressed on most normal epithelial cells and gastrointestinal carcinomas and functions as a homotypic calcium-independent cell adhesion molecule. The antigen is being used as a target for immunotherapy treatment of human carcinomas. Mutations in this gene result in congenital tufting enteropathy.
基因ID	4072
基因名	EPCAM
Swiss	P16422 (https://www.uniprot.org/uniprotkb/P16422/entry)
别名	EpCAM (YD18511),EpCAM (YD18511) Rabbit mAb,EPCAM,Adenocarcinoma-associated antigen,Cell surface glycoprotein Trop-1,Epithelial cell surface antigen,Epithelial glycoprotein,Epithelial glycoprotein 314,KSA 1/4 antigen,KSA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)