

Aurora A (YD18466) Rabbit mAb

货号: **AYD14685**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	46kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HT-29
细胞定位	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, spindle pole, centriole, Cell projection, neuron projection, cilium, cilium basal body, Basolateral cell membrane
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a cell cycle-regulated kinase that appears to be involved in microtubule formation and/or stabilization at the spindle pole during chromosome segregation. The encoded protein is found at the centrosome in interphase cells and at the spindle poles in mitosis. This gene may play a role in tumor development and progression. A processed pseudogene of this gene has been found on chromosome 1, and an unprocessed pseudogene has been found on chromosome 10. Multiple transcript variants encoding the same protein have been found for this gene.
基因ID	6790
基因名	AURKA

Swiss	O14965
别名	Aurora A (YD18466)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn