

— ABLYBIO, Help Your Research



# IL-6 (YD18277) Rabbit mAb

货号: **AYD14717**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	24kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	Cytokine with a wide variety of biological functions in immunity, tissue regeneration, and metabolism (Probable). Binds to IL6R, then the complex associates to the signaling subunit IL6ST/gp130 to trigger the intracellular IL6-signaling pathway (PubMed:8910279). The interaction with the membrane-bound IL6R and IL6ST stimulates 'classic signaling', whereas the binding of IL6 and soluble IL6R to IL6ST stimulates 'trans-signalng'. Alternatively, 'cluster signaling' occurs when membrane-bound IL6:IL6R complexes on transmitter cells activate IL6ST receptors on neighboring receiver cells (PubMed:27893700)
基因ID	16193
基因名	IL6
Swiss	P08505 ( <a href="https://www.uniprot.org/uniprotkb/P08505/entry">https://www.uniprot.org/uniprotkb/P08505/entry</a> )
别名	IL-6 (YD18277),IL-6 (YD18277) Rabbit mAb,IL6,B-cell hybridoma growth factor,Interleukin HP-1,IL-6

## 产品验证

Recombinant Mouse IL-6 protein



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)