

PDCD6/ALG-2 (YD16835) Rabbit mAb

货号: **AYD14779**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	22kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Endoplasmic reticulum membrane, Cytoplasmic vesicle, COPII-coated vesicle membrane, Cytoplasm, Nucl eus, Endosome
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a calcium-binding protein belonging to the penta-EF-hand protein family. Calcium binding is important for homodimerization and for conformational changes required for binding to other protein partners. This gene product participates in T cell receptor-, Fas-, and glucocorticoid-induced programmed cell death. In mice deficient for this gene product, however, apoptosis was not blocked suggesting this gene product is functionally redundant. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is also located on the short arm of chromosome 5.
基因ID	10016

基因名	PDCD6
Swiss	O75340
别名	PDCD6/ALG-2 (YD16835)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn