

— ABLYBIO, Help Your Research



G3BP (YD19076) Rabbit mAb

货号: **AYD14815**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	52kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Cytoplasm, cytosol, Perikaryon, Stress granule, Nucleus
纯化	亲和纯化




抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes one of the DNA-unwinding enzymes which prefers partially unwound 3'-tailed substrates and can also unwind partial RNA/DNA and RNA/RNA duplexes in an ATP-dependent fashion. This enzyme is a member of the heterogeneous nuclear RNA-binding proteins and is also an element of the Ras signal transduction pathway. It binds specifically to the Ras-GTPase-activating protein by associating with its SH3 domain. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined.
基因ID	10146
基因名	G3BP1
Swiss	Q13283 (https://www.uniprot.org/uniprotkb/Q13283/entry)
别名	G3BP (YD19076),G3BP (YD19076) Rabbit mAb,G3BP1,ATP-dependent DNA helicase VIII,GAP SH3 domain-binding protein 1,G3BP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)