

— ABLYBIO, Help Your Research



Interferon beta (YD16512) Rabbit mAb

货号: **AYD14849**

产品信息

| | |
|-------|---|
| 反应 | Mouse |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IP |
| 推荐浓度 | |
| 理论分子量 | 22kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Recombinant Mouse IFN-beta Protein |
| 细胞定位 | Secreted |
| 纯化 | 亲和纯化 |

抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|---|
| 研究背景 | Enables cytokine activity and type I interferon receptor binding activity. Involved in B cell proliferation; cellular response to virus; and type I interferon signaling pathway. Acts upstream of or within several processes, including cellular response to organic cyclic compound; defense response to other organism; and negative regulation of osteoclast differentiation. Located in extracellular space. Is expressed in embryo and liver. Human ortholog(s) of this gene implicated in COVID-19. Orthologous to human IFNB1 (interferon beta 1). |
| 基因ID | 15977 |
| 基因名 | ifnb1 |
| Swiss | P01575 (https://www.uniprot.org/uniprotkb/P01575/entry) |
| 别名 | Interferon beta (YD16512),Interferon beta (YD16512) Rabbit mAb,ifnb1,ifb,ifnb |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)