

— ABLYBIO, Help Your Research



# GATA4 (YD16510) Rabbit mAb

货号: **AYD14850**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	45kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,Mouse testis,Rat stomach
细胞定位	Nucleus
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the GATA family of zinc-finger transcription factors. Members of this family recognize the GATA motif which is present in the promoters of many genes. This protein is thought to regulate genes involved in embryogenesis and in myocardial differentiation and function, and is necessary for normal testicular development. Mutations in this gene have been associated with cardiac septal defects. Additionally, alterations in gene expression have been associated with several cancer types. Alternative splicing results in multiple transcript variants.
基因ID	2626
基因名	GATA4
Swiss	P43694 ( <a href="https://www.uniprot.org/uniprotkb/P43694/entry">https://www.uniprot.org/uniprotkb/P43694/entry</a> )
别名	GATA4 (YD16510),GATA4 (YD16510) Rabbit mAb,GATA4,GATA-binding factor 4

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)