

— ABLYBIO, Help Your Research



FOXG1 (YD16286) Rabbit mAb

货号: **AYD14873**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	52kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,SH-SY5Y,Mouse brain,Mouse spleen,Rat brain
细胞定位	Nucleus
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This locus encodes a member of the forked-head transcription factor family. The encoded protein, which functions as a repressor, may play a role in brain development. Mutations at this locus have been associated with Rett syndrome.
------	---

基因ID	2290
基因名	FOXG1
Swiss	P55316 (https://www.uniprot.org/uniprotkb/P55316/entry)
别名	FOXG1 (YD16286),FOXG1 (YD16286) Rabbit mAb,FOXG1,Brain factor 1,Brain factor 2,Forkhead box protein G1A,Forkhead box protein G1B,Forkhead box protein G1C,Forkhead-related protein FKHL1,Forkhead-related protein FKHL2,Forkhead-related protein FKHL3,FKH2,FKHL1,FKHL2,FKHL3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)