

— ABLYBIO, Help Your Research



Chk1 (YD12192) Rabbit mAb

货号: **AYD14895**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC IP
推荐浓度	
理论分子量	54kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,HepG2,C6
细胞定位	Nucleus, Chromosome, Cytoplasm, cytoskeleton, microtubule organizing center, centrosome
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene belongs to the Ser/Thr protein kinase family. It is required for checkpoint mediated cell cycle arrest in response to DNA damage or the presence of unreplicated DNA. This protein acts to integrate signals from ATM and ATR, two cell cycle proteins involved in DNA damage responses, that also associate with chromatin in meiotic prophase I. Phosphorylation of CDC25A protein phosphatase by this protein is required for cells to delay cell cycle progression in response to double-strand DNA breaks. Several alternatively spliced transcript variants have been found for this gene.
基因ID	1111
基因名	CHEK1
Swiss	O14757 (https://www.uniprot.org/uniprotkb/O14757/entry)
别名	Chk1 (YD12192), Chk1 (YD12192) Rabbit mAb, CHEK1, CHK1 checkpoint homolog, Cell cycle checkpoint kinase, Checkpoint kinase-1, CHK1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)