

— ABLYBIO, Help Your Research



IL-1 alpha (YD16122) Rabbit mAb

货号: **AYD14909**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	31kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Nucleus, Cytoplasm, Secreted
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a member of the interleukin 1 cytokine family. This cytokine is a pleio tropic cytokine involved in various immune responses, inflammatory processes, and hematopoiesis. This cytokine is produced by monocytes and macrophages as a proprotein, which is proteolytically processed and released in response to cell injury, and thus induces apoptosis. This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. It has been suggested that the polymorph ism of these genes is associated with rheumatoid arthritis and Alzheimer's disease.
基因ID	3552
基因名	IL1A
Swiss	P01583 (https://www.uniprot.org/uniprotkb/P01583/entry)
别名	IL-1 alpha (YD16122),IL-1 alpha (YD16122) Rabbit mAb,IL1A,Hematopoietin-1,IL1F1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)